
**Komplet vsebuje**

- a) Priključek 2 x
- b) Zbiralnik nosnega izcedka 2 x
- c) Povezovalna cev 1 x
- d) Filter 5 x

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
**Accessory kit**


- a) Connector 2 x
- b) Plastic mucus container 2 x
- c) Flexible tubing 1 x
- d) Replacing filter 5 x

**Why is it advisable to use a nasal aspirator for babies?**


Infants up to one year of age only breathe through nose and they are not able to blow out their own mucus. Eliminated mucus in the baby nose makes difficulties at breathing. The nasal aspirators ME126 helps at removing the mucus out from infant nozzle. Clean and transitional nostrils allow an infant a peaceful sleep and easier breast-feeding.

**Nasal aspirator ME126 assembling**


1. Insert the connector into the tube

2. Insert into connector clean filter

3. Connector and filter cover with container

**Cleaning of the nasal aspirator**

- After each use of the nasal aspirator ME126, disassemble and wash it with hot water and soap.
- After use always store already cleaned nasal aspirator.
- If the nasal aspirator is used continuously we recommend that connectors and plastic containers are boiled in the water.

**1. Classical aspiration through the mouth**

- a) Wash your hands.
- b) Assemble the nasal aspirator and insert a new clean filter.
- c) Place the baby on the table and wedges his head.
- d) Place free connector in your mouth.
- e) Top of the plastic mucus container gently insert into the nostril.
- f) With gentle breath suck baby mucus in to plastic mucus container.
- g) The procedure should be repeated several times if necessary.
- h) The procedure repeat with the other nostril.

Prior inserting saline solution into the nostrils facilitates secretion and absorption of mucus, so the purchase of nebulizer ME100 is recommended.

[www.medikoel.com](http://www.medikoel.com)

**2. Aspiration with help of piston nebulizer ME100**

- a) Wash your hands.
- b) Assemble the nasal aspirator and insert a new clean filter.
- c) Place the baby on the table and wedges his head.
- d) From the nebulizer ME100 remove filter cap.
- e) Connect the free connector on the inlet of the nebulizer ME100.
- f) Switch on the nebulizer ME100 by pressing the switch in position '1'.
- g) Top of the plastic mucus container gently insert into the nostril..
- h) If the mucus container is full replace it with the clean one.
- i) The procedure repeat with the other nostril.
- j) After the aspiration is finished stop nebulizer ME100 by pressing the switch in position '0'.

Prior inserting saline solution into the nostrils facilitates secretion and absorption of mucus, therefore the prior use of nebulizer ME100 is recommended.