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INSTRUCTIONS FOR USE

PLEASE REFER TO THE INSTRUCTIONS FOR USE

THESE INSTRUCTIONS ARE IMPORTANT. PLEASE KEEP THEM FOR FUTURE REFERENCE.

Dear customer, Thank you for choosing Pic Solution's LiteRAPID, an automatic blood pressure monitor designed and manufactured using the most up-to-date technology. Before using the device, we highly recommend that you read the short list of warnings contained in this manual, to make sure you have correctly understood how to use it.

INTRODUCTION

The regular use of a Blood Pressure Monitor allows you to monitor blood pressure continuously, but it should not take the place of regular check-ups by your doctor. We recommend that you visit your doctor regularly for general check-ups and for more detailed information on blood pressure. Blood pressure measurements determined with this unit are equivalent to those obtained by a trained observer using cuff/stethoscope auscultation method, within the limits prescribed by the international ESH Validation protocol.


WARRANTY

The product is guaranteed for 3 years against all conformity defects under the normal utilisation conditions provided in the user instructions. The warranty shall not therefore apply in the case of damages caused by improper use, wear or accidental events. Batteries, cuff and wearing parts are not included in the warranty. In case of malfunction, contact exclusively the authorised Pic Technical Customer Service.

RAPID-TECH Comfort Experience™

This unit uses the oscillometric method to detect your blood pressure. RAPID-TECH Comfort Experience™ is the latest generation technology chosen by PIC for the blood pressure monitors of the line RAPID. It allows fast and, therefore, less annoying monitoring because the measurement is performed already during the cuff inflation phase. Pic LiteRAPID has been clinically validated in accordance with the ESH (European Society of Hypertension) protocol, as a guarantee of its accuracy and precision.

NAME/FUNCTION OF EACH PART

List of included items: Monitor, M-L size Cuff (REF 02007827000000), 4xAAA batteries, User manuals, Carrying bag.

ABOUT BLOOD PRESSURE

What is blood pressure?

Blood pressure is the pressure exerted on the artery tube while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

Why measure your blood pressure?


Among the various health problems afflicting modern people, problems associated with high blood pressure are by far the most common. High blood pressure's dangerously strong correlation with cardiovascular diseases and high morbidity has made measuring blood pressures a necessity of identifying those at risk.

Blood Pressure Fluctuation

Blood pressure fluctuates all the time! You should not be over worried if you encountered two or three measurements at high levels. Blood pressure changes over the month and even throughout the day. It is also influenced by the season, temperature and activities.

Classification of high blood pressure

The out-of-office blood pressure classification published by ESH/ESC in 2013 is as follows:

Note: It is not intended to provide a basis of any type of rush toward emergency conditions/diagnosis based on the scheme; the scheme is meant only to discriminate between the different levels of blood pressure.

WARNING

Only a physician can tell your normal BP range. Please contact a physician if the average of your measuring results falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

Irregular Heartbeat Detector

This blood pressure monitor is equipped with an algorithm for the detection of irregular heartbeats (IHB) and is capable of reporting irregularities when the heart rhythm changes while the monitor is measuring systolic and diastolic blood pressure. During each measurement, the monitor records all the pulse intervals and calculate the average; if there are two or more pulse intervals, the difference between each interval and the average is more than the average value of $\pm 25\%$, or there are four or more pulse intervals, the difference between each interval and the average is more than the average value of $\pm 15\%$, the irregular heartbeat symbol appears on the display when the measurement results are appear.

WARNING

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. We wish to underline that the device does not replace a visit to the cardiologist, but it is useful to detect any irregularities in advance.

Why the blood pressure I get from the hospital is different from home?


The blood pressure is different even during 24 hour because of the weather, emotion, exercise etc, specially the "white coat" in hospital which makes the results are higher than the ones at home.

Is the result the same if measuring on the right arm?

Measurements can be taken on either arm; results can however differ sometimes so it is always advisable to use the same arm.

LCD display explanation

At starting check the completeness of the LCD segments

SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic pressure	Maximum/systolic pressure result
DIA	Diastolic pressure	Minimum/diastolic pressure result
▼	Pulse	Pulse/minute
♥	Deflating	The arm cuff is deflating
88	Memory	If "M" shows, the displayed measurement values is from the memory
mmHg	mmHg	Blood pressure unit of measurement mmHg
⦿ + ⦿	Low battery	Batteries are low and need to be replaced
♥	Irregular heartbeat	Irregular heartbeat

👤 👤	User 1/User 2	Select User 1/User 2
📊	Grade	Blood pressure level indicator
M 08:08 08:08	Date	"MM" shows the month, "DD" shows the day
Avg	Average	The average of last 3 measurements taken within 30 minutes.

The Choice of Power Supply

Dry cell powered mode: 4 dry cells 1.5V AAA; 6V DC.

WARNING

To protect the device and your safety and to achieve optimal performance, use only the dry cells and adaptors specified in this manual (4 x 1.5V AAA dry cells; 6V DC).

Installing and Replacing the Dry cells

- Slide off the battery cover.
- Install the batteries by matching the correct polarity, as shown by the marks inside the battery compartment.
- Replace the cover.

Replace the batteries whenever the below happen

Is displayed **⦿ + ⦿** the display dims, the display does not light up

WARNING

- If the device is not used for a while, please remove the batteries.
- The dry cells used are harmful for the environment. Dispose of in accordance with the instructions in this manual (EU Directive 2006/66/EC).
- Do not dispose of batteries in fire. Batteries may explode or leak.
- For instructions on battery disposal please refer to the indication reported in this user manual.
- There are no user serviceable parts inside. Dry cells or damage from old dry cells is not covered by warranty.
- Only use brand batteries. Always replace with new dry cells together. Use dry cells of the same brand and same type. Other types of dry cells can cause explosions or fires.
- The battery can explode if handled carelessly. Never recharge, remove or position the dry cells near sources of heat. Do not incinerate.
- Rechargeable batteries are not recommended.


Setting Date and Time

It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each measurement that is stored in the memory. (year: 2020-2060, time: 12H/24H). In the setting mode, use the "MEM" button to change the value and the "SET" button to confirm your choice.

- When the unit is off, press "SET" button, then hold pressing "SET" for 3 seconds to enter the mode for year setting.
- Press the "MEM" to change the [YEAR].
- When you get the right year, press "SET" to set down and turn to next step.
- Repeat step 2 and 3 to set the [MONTH] and [DAY].
- Then it will enter time setting mode. There are two kinds of time format "12" hours and "24" hours. Press "MEM" button to change the time format. Press "SET" button to confirm the time format.
- Repeat step 2 and 3 to set the [HOUR] and [MINUTE].
- After confirming the time, the LCD will display "00E" and the monitor will shut off.

Wear the cuff

- Remove any watches, jewellery, etc. from the left arm. Note: if a doctor has diagnosed limited circulation in the left arm, use the right arm.
- Roll up the sleeve to expose the skin. Be sure that the sleeve is not tight.
- Plug in the cuff connecting tube into the unit.
- Wrap the arm cuff, leaving the end of the cuff around the cuff.

Helpful Tips

- Here are a few helpful tips to help you obtain more accurate readings:
- Do not measure your blood pressure immediately after consuming a large meal.
 - To obtain more accurate readings, please wait one hour after meal before measuring.
 - Do not smoke or drink tea, coffee or alcohol before measuring your blood pressure.
 - Do not measure your blood pressure immediately after taking a bath; wait at least 20 minutes.
 - You should not be physically tired or exhausted while taking a measurement.
 - Do not measure your blood pressure when you want to discharge urine.
 - It is important that you relax during measurement.
 - Do not take measurements if you are under stress or tense.
 - Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.
 - Rest for 5 minutes before measuring.
 - Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
 - Take the measurement in a silent room.
 - Do not move and do not talk when taking a measurement
 - Do not cross your legs and keep your feet on the ground.
 - Keep your back against the backrest of the chair.
 - For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

Select the User

When the monitor is off, press and hold the MEM button to enter user setting mode, the user ID will blink. Press MEM button again, select the user ID and press SET button to save; then "user ID + done" will be shown and the device will turn off.

Start the measurement

- When the monitor is off, press the **START/STOP** **▶** to turn on the monitor and start the measurement.

If the **🔊** icon appears, please, read carefully the Helpful Tips section and repeat the measurement. **Notes:** To stop the measurement at any time, press the START/STOP button. The cuff will quickly deflate. Store the cuff carefully.

Recall the measurements

- When the Blood pressure monitor is off, press the "MEM" button to see the average of the last 3 recordings taken within 30 minutes. If the recorded measurements are less than 3 as well as the last 3 recordings are taken in more than 30 minutes, you will see the latest reading.
- Press "MEM" or "SET" to see the recording wanted.

WARNING

The most recent recording (1) is shown first. Every time a new recording is made, the first (1) record is assigned to it. All the others are moved by one figure (e.g. 2 becomes 3, and so on), whereas the last record (60) is deleted from the list.

Delete the measurements

If you did not get the correct measurement, you can delete all results by following below steps.

- Hold pressing "MEM" for 3 seconds when the monitor is in the memory recall mode. The display will read "del All" and the User number. Press "SET" to confirm deleting and the monitor will turn off.
- If you don't want to delete the records, press "START/STOP" to escape.
- If there is no record the display will show "-,-".

Maintenance

- In order to get the best performance from your monitor, please follow the below instructions.
- Put in a dry place and avoid the sunshine.
 - Avoid touching water, clean it with a dry cloth in case.
 - Avoid the intense shaking and collision.
 - Avoid the dusty and unstable temperature environment.
 - Using the wet clothing to remove the dirt.
 - Avoid washing the cuff.

Calibration and service

The accuracy of this monitor has been carefully tested and is has been designed for a long life. Indeed the monitor can main-

tain the safety and performance characteristics for a minimum of 10,000 measurements or five years of normal use. It is generally recommended to check the unit every two years to ensure correct functioning and accuracy. Please consult the Contact information in this user manual to get in touch with your local authorized Service Center.

WARNINGS

- The patient and the operator are the same person.
- The sphygmomanometer has been clinically tested as required by ISO 81060-2.
- This device is intended for non-invasive testing and monitoring of arterial blood pressure in adults at home. This device is not suitable for use on newborns, children or pregnant women. The device is not validated to be used on the above listed patient and it could be dangerous. The use on teenagers shall be evaluated by physician.
- Use of the device is contraindicated in women who may be pregnant or are pregnant or who suffer from pre-eclampsia. Besides providing inaccurate readings, the effects of this device on the fetus are unknown.
- The device is not suitable for use if your arm circumference is not included in the range 22cm-42cm (8 1/2" - 16 1/2").
- People with reduced cognitive abilities shall use the device under supervision.
- Keep out of reach of children, people with reduced cognitive abilities or pets: the equipment contains small parts which can cause suffocation or internal lesions if swallowed. The tubes of the arm cuff may pose a strangulation hazard.
- Blood pressure measurement is not suitable in cases of serious arteriosclerosis (hardening of the arteries).
- If the arm cuff pressure reaches 300 mmHg it deflate automatically. If the arm cuff does not deflate when pressure reaches 300 mmHg or if you feel discomfort during the measurement take it off your arm and press the START/STOP button to stop inflation.
- Do not confuse self-monitoring with self-diagnosis. By using this device you can monitor your blood pressure, but it does not replace your doctor. Begin or end medical treatment based solely on medical advice. If you are taking medication, consult your physician to determine the most appropriate time to test your blood pressure. Never change a prescribed medication without consulting your Physician.
- If an Irregular Heartbeat (IHB) caused by a common arrhythmia is detected while taking a blood pressure reading, a signal will be displayed. Under this condition, the electronic sphygmomanometer can keep functioning, but the results may not be accurate. It is suggested that you consult with your physician for accurate assessment. When the device is used for measurements on patients with common arrhythmias, such as atrial or ventricular premature beats or atrial fibrillation, a deviation may occur. Please consult your doctor about the result.
- Please consult your physician before using the device if you have active implanted medical devices.
- This device cannot be used with HF surgical equipment at the same time.
- This device is not intended to be used during patient transportation outside a healthcare facility.
- It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- The user must check that the equipment functions safely and see that it is in proper working condition before being used.
- Manufacturer will make available on request circuit diagrams, component parts list etc. only to qualified person.
- Too frequent measurements can cause injury to the PATIENT due to blood flow interference. Please check by observation of the limb concerned) that operation of LiteRAPID does not result in prolonged impairment of the circulation of the blood of the PATIENT.
- This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's arm and fingers will become anaesthetic, swollen and even purple due to a lack of blood.
- Check that the device tubes are not kinked while using the device: kinked tubes can generate lack of blood in the extremities. The equipment includes delicate components and must be handled with care.
- Please use the device under the environment which was provided in the instruction manual. If this is not complied with, correct operation and life of the device may be compromised.
- Do not use if allergic to polyester, nylon or plastic.
- During using, the patient will contact with the cuff. The materials of the arm cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.
- The application of the CUFF over a wound can cause further injury.
- The application of the CUFF and its pressurization on any limb where intravenous access or therapy, or an arterio-venous (A-V) shunt, is present because of temporary interference to blood flow and could result in injury to the PATIENT
- Applying the ARM CUFF and inflating it on an arm with a mastectomy on its side may result in injury, in addition, the testing may not be accurate.
- Do not inflate the cuff on the same limb which other monitoring ME equipment is applied around simultaneously, because this could cause temporary loss of function of the simultaneously used monitoring ME equipment.
- The device is a non-AP/APG device and not suitable for use in the presence of anaesthetic mixtures that are inflammable with air, oxygen or nitrous oxide.
- The operator must not touch the battery terminals and the patient at the same time.
- Please use ACCESSORIES and detachable parts specified/authorised by the MANUFACTURER. Otherwise, it may cause damage to the unit or danger to the user/patient.
- We recommend checking performance every two years or after repairs. The device must be checked again in case of damage due to impact (such as falls) or exposure to liquids and/or extreme temperatures (hot/cold) or extreme humidity variations.
- For public place use, please contact your local distributor or the Italian toll-free number 800 900 080 to get support how to proceed.
- Please dispose of accessories, detachable parts, and the ME EQUIPMENT according to the local guidelines.
- In case of problems concerning performances, settings, maintenance or use of the device, contact your local distributor or the Italian toll-free number 800 900 080. **Warning!** Do not service or repair the appliance while it is being used! Any maintenance/service activity/operation must be done solely by a PIC Service Center. No modification of the equipment are allowed.
- When not in use, store the equipment in a dry area and guard against dust, heat, humidity and direct sunlight. Do not place heavy objects on the travel bag.
- If the device is stored in an environment at the minimum/maximum storage temperature, wait 30 minutes for the device to heat up/cool down to the right level of temperature for use.

This section includes a list of error messages and frequently asked questions for problems you may encounter with your sphygmomanometer. If the products not operating as you think it should, check here before arranging for servicing.

PROBLEM	SYMPTOM	CHECK THIS	REMEDY
No power	Display will not light up.	Batteries are exhausted.	Replace with new batteries
		Batteries are inserted incorrectly.	Insert the batteries correctly
Batteries exhausted	The display dims or reads ⦿ + ⦿	Batteries are exhausted.	Replace with new batteries
	E 1 shows	The cuff is not secure or too tight	Refasten the cuff and then measure again.
	E 2 shows	The monitor detected motion while measuring	Refasten the cuff and then measure again.
	E 3 shows	the measurement process does not detect the pulse signal.	Relax for a moment and then measure again.
Error message	E 4 shows	The treatment of the measurement failed.	Relax for a moment and then measure again.
	EEx shows on the display.	A calibration error occurred.	Retake the measurement. If the problem persists, contact the dealer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions.
Error message	OUT	Out of range measurement.	Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your Physician.

SPECIFICATIONS

- Power supply:** Dry cell powered mode: 4 dry cells 1.5V AAA; 6V DC (~105 tests).
- Display:** LCD, VA, 52 x 58 mm
- Measurement mode:** Oscillometric inflation -type measurement mode
- Measurement range:** Rated arm cuff pressure: from 0mmHg to 299mmHg
Measurement pressure: SYS: from 60 to 230mmHg, DIA: from 40 to 130mmHg
pulse value: from 40 to 199 beats/minute
Pressure: (from 5°C to 40°C) \pm 3mmHg
pulse value: \pm 5%
- Accuracy:** Temperature: from 5°C to 40°C. Relative humidity from 15% to 90%
Atmospheric pressure: 70kPa to 106kPa
- Normal working condition:** Temperature: from -20°C to 60°C; Relative humidity \leq 93%
Atmospheric pressure: 50kPa to 106kPa
- Storage and transport conditions:** Temperature: from -20°C to 60°C; Relative humidity \leq 93%
Atmospheric pressure: 50kPa to 106kPa
- Arm circumference:** About 22cm-42cm
- Net Weight:** Approx. 200g (excluding the dry cells)
- External dimensions:** Approx. 102mm x 107mm x 40mm
- Accessories:** REF 02007827000000
- Mode of Operation:** Continuous operation
- Useful life:** 10,000 Measurements
- Degree of protection:** Type BF applied part (cuff)
- Protection against ingress of water:** IP21 degree of protected against solid foreign objects of 12.5mm ϕ and greater and protected against vertically falling water drops A01
- Software version:** A01

Contact Information

For further information about our products, please visit our website www.picsolution.com where you can receive assistance and address the most common issues.

ELECTROMAGNETIC COMPATIBILITY

Electro-medical devices necessitate specific precautions in terms of electromagnetic compatibility and must be installed and prepared for use according to the following EMC Information. Electromagnetic radio frequency interference (EMC/RFI): the readings may be subject to interference if the device is used in an electromagnetic field with a radio frequency of \leq 10 Volts/metre but the instrument's effectiveness will not be compromised permanently. Keep the thermometer at a distance of at least 1 m from RF transmitters in order to avoid radio frequency interference.

Manufacturer's declaration and guidelines - electromagnetic emissions		
Phenomenon	Professional healthcare facilities ^{*)}	HOUSEHOLD ENVIRONMENT ^{*)}
Conducted and radiated RF EMISSIONS	^{*)}	CISPR 11 Group 1 Class B

Harmonic distortion and voltage fluctuations/flicker	Not Applicable
--	----------------

a) The device was designed to be used in a domestic setting or in professional healthcare facilities and can be used exclusively in hospital wards and rooms allocated for outpatient treatment/diagnosis in hospitals and clinics. The more restrictive acceptance limits of Group 1 Class B (CISPR 11) have been considered and applied.

Manufacturer's declaration and guidelines - Electromagnetic immunity - Port: Housing			
Phenomenon		Immunity test levels	
		Professional healthcare facilities	HOUSEHOLD ENVIRONMENT ^{*)}
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	\pm 8kV contact \pm 2kV, \pm 4kV, \pm 15kV air	
Radiated RF EM fields	IEC 61000-4-3	^{*)}	10 V/m ^{*)} 80MHz - 2.7GHz 80% AM at 1kHz
Magnetic fields at NOMINAL grid frequency	IEC 61000-4-8	30 A/m ^{*)} 50Hz or 60Hz	

^{*)} The device was designed to be used in a domestic setting or in professional healthcare facilities, exclusively in hospital wards and rooms allocated for outpatient treatment/diagnosis in hospitals and clinics. The most restrictive acceptable limits were considered and applied in the immunity tests.

^{*)} Before the application of modulation.

^{*)} This test level presupposes a minimum distance of at least 15 cm between the LiteRAPID REF 02022533000000 and sources of magnetic fields at grid frequency.

	"PLEASE REFER TO THE INSTRUCTIONS FOR USE" symbol		"TYPE BF APPLIED PART" symbol
	Symbol for "COMPLIES WITH MDD 93/42/EEC as amended by directive: 2007/47/EC"		Storage and transport humidity symbol
	"MANUFACTURER" symbol		Storage and transport temperature symbol
	Symbol for "SERIAL NUMBER"		Warning
	Symbol for "DIRECT CURRENT"		Lot number
	Latex free		Catalogue number
	Protected against solid foreign objects of 12.5mm ϕ and greater and protected against vertically falling water drops.		Atmospheric pressure

Appliance date of production: the second pair of digits in the lot number identifies the year, the third pair identifies the month. (LOT) 1016100001; 16 = 2016; 10 = October)

THE PRODUCT COMPLIES WITH:

EN 1060-3: 1997/A2:2009 Non-invasive sphygmomanometers Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems.
EN 1060-4:2004 Non-invasive sphygmomanometers Part 4: Procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.
IEC 80601-2-30:2009 + A1:2013 Medical electrical equipment - Part 2: Special regulations concerning the fundamental safety and essential performance of non-invasive automatic sphygmomanometers.

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disposed device is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

THIS PRODUCT COMPLIES WITH EU DIRECTIVE 2006/66/EC

The crossed bin symbol on the batteries indicates that, at the end of their life, they must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for bringing batteries, at the end of their life, to the appropriate collection facilities in order to facilitate treatment and recycling. The correct separate waste disposal of your old toy will help prevent potential negative consequences for the environment and human health and favour the recycling of the materials of which the product is made. Illegal disposal of the product by the user will damage the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

Pouch cleaning symbols

- Wash by hand in cold water,
- Do not bleach,
- Do not tumble dry,
- Line drying,
- Do not iron,
- Do not dry clean.

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