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INSTRUCTIONS FOR USE

PLEASE REFER TO THE INSTRUCTIONS FOR USE

THESE INSTRUCTIONS ARE IMPORTANT. PLEASE KEEP THEM FOR FUTURE REFERENCE.

Dear customer,
Thank you for choosing Pic Solutions maxiRAPID, an automatic blood pressure monitor designed and manufactured according to the most up-to-date technologies. Before using the device, we highly recommend that you read the short list of warnings contained in this manual, to make sure you have correctly understood how to use it.

INTRODUCTION


The regular use of a Blood Pressure Monitor allows you to monitor blood pressure continuously, but it should not take the place of regular check-ups by your doctor. We recommend that you visit your doctor regularly for general check-ups and for more detailed information on blood pressure. Blood pressure measurements determined with this unit are equivalent to those obtained by a trained observer using cuff/ stethoscope auscultation method, within the limits prescribed by the International Protocol (ESH Validation).

This unit uses the oscillometric method to detect your blood pressure. The measurement is performed already during the cuff inflation phase. Pic maxiRAPID has been clinically validated in accordance with the ESH (European Society of Hypertension) protocol, as a guarantee of its accuracy and precision.

WARRANTY

The product is guaranteed for 3 years against all conformity defects under the normal utilisation conditions provided in the user instructions. Therefore, the warranty shall not apply in case of damages caused by improper use, wear, or accidental events. Batteries, cuff and wearing parts are not included in the warranty. In case of malfunction, contact exclusively the authorised Pic Technical Support Centre.

NAME/FUNCTION OF EACH PART

List of included items: Blood pressure monitor, M-L size cuff (REF 0X007827000000 X= 0,2), 4xAAA batteries, user manual, carrying bag. Features:
• 61x93 mm LCD display,
• Two memory areas
• 100 recorded measurements per area,
• Inflation-type measurement.

ABOUT BLOOD PRESSURE

What is blood pressure?
Blood pressure is the pressure exerted on the artery tube while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

Why measure your blood pressure?


Among the various health problems afflicting modern people, problems associated with high blood pressure are by far the most common. High blood pressure's dangerously strong correlation with cardiovascular diseases and high morbidity has made measuring blood pressures a necessity of identifying those at risk.

Blood Pressure Fluctuation

Blood pressure fluctuates all the time! You should not be over worried if you encountered two or three measurements at high levels. Blood pressure changes over the month and even throughout the day. It is also influenced by the season, temperature and activities.

Classification of blood pressure

The out-of-office blood pressure classification published by ESH/ESC in 2013 is as follows:

Note: It is not intended to provide a basis of any type of rush toward emergency conditions/diagnosis based on the scheme: the scheme is meant only to discriminate between the different levels of blood pressure.

WARNING
Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

Irregular Heartbeat Detector

This blood pressure monitor is equipped with an algorithm for Irregular Heartbeat (IHB) Detection. An irregular heartbeat is detected when a heartbeat rhythm varies while the unit is measuring the systolic and diastolic blood pressure. During each measurement, the monitor records all the pulse intervals and calculate the average; if there are two or more pulse intervals, the difference between each interval and the average is more than the average value of ±25%, or there are four or more pulse intervals, the difference between each interval and the average is more than the average value of ±15%, then the irregular heartbeat symbol will appear on the display with the measurement result.


WARNING
The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. We wish to underline that the device does not replace a visit to the cardiologist, but it is useful to detect any irregularities in advance.

Why the blood pressure I get from the hospital is different from home?
The blood pressure is different even during 24 hour because of the weather, emotion, exercise etc, specially the "white coat" in hospital which makes the results are higher than the ones at home.

Is the result the same if measuring on the right arm?
Measurements can be taken on either arm; results can however differ sometimes so it is always advisable to use the same arm.

LCD display explanation

At starting check the completeness of the LCD segments

SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic pressure	Maximum/systolic pressure result
DIA	Diastolic pressure	Minimum/diastolic pressure result
Pul/min	Pulse	Pulse/minute
▼	Deflating	The arm cuff is deflating
MEM	Memory	If "M" shows, the displayed measurement values is from the memory
mmHg	mmHg	Blood pressure unit of measurement mmHg
Lo +	Low battery	Batteries are low and need to be replaced
♥	Irregular heartbeat	Irregular heartbeat
I	Grade	Blood pressure level indicator
00:00	Current Time	Year/Month/Day, Hour/Minute
AVG	Average	The average of last 3 measurements taken within 30 minutes.
A	User A	Start measurement for User A
B	User B	Start measurement for User B
♥	Heartbeat	Heartbeat detection during the measurement

The Choice of Power Supply
1. Battery powered mode: 4 batteries 1.5V AAA; 6V ≈ 1A
2. Power supply via mains adapter: Input: 100 - 240V~, 50 / 60Hz, 0.2A max, Output: 6V ≈ 1A
Can be powered only with the "Adapter for blood pressure monitor" (REF 02000520000000 (BLJ06L060100P-V)) with European plug and the "Adapter for blood pressure monitor UK" (REF 02010310000000 (BLJ06L060100P-B)) with UK plug (Not included). Please unplug the adaptor to depart from the using utility power.
NB: The adaptor location is located on the right side of the monitor. Place NO obstacles on the right side for easy pull-off adaptor.

WARNING
In order to protect your device and yourself and to get the best performances, use only the batteries and adaptors according to the information provided in this manual.


Installing and Replacing the Batteries
1. Slide off the battery cover.
2. Install the batteries by matching the correct polarity, as shown by the marks inside the battery compartment.
3. Replace the cover.

Replace the batteries whenever the below happen
Is displayed **Lo +**, the display dims, the display does not light up.

WARNING
• If the device is not used for a while, please remove the batteries.
• The old battery is harmful to the environment, so please disposal accordingly to the indication reported in this user manual (EU Directive 2006/66/EC).
• Do not dispose of batteries in fire. Batteries may explode or leak.
• For instructions on battery disposal please refer to instructions.
• There are no user serviceable parts inside. Batteries or damage from old batteries is not covered by warranty.
• Only use brand batteries. Always play with new batteries together. Use batteries of the same brand and same type. Other types of batteries can cause explosions or fires.
• The battery can explode if handled carelessly. Never recharge, remove or position the batteries near sources of heat. Do not incinerate.
• Rechargeable batteries are not recommended.


Setting Date and Time
It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each measurement that is stored in the memory. (year: 2015-2050, time: 12H/24H). In the setting mode, use the "MEM" button to change the value and the "SET" button to confirm your choice.
1. When the unit is off, hold pressing "SET" for 3 seconds to enter the mode for year setting.
2. Press the "MEM" to change the [YEAR].
3. When you get the right year, press "SET" to set down and turn to next step.
4. Repeat step 2 and 3 to set the [MONTH] and [DAY].
5. Then it will enter time setting mode. There are two kinds of time format "12" hours and "24" hours. Press "MEM" button to change the time format. Press "SET" button to confirm the time format.
6. Repeat step 2 and 3 to set the [HOUR] and [MINUTE].
7. After confirming the time, the LCD will display "OnE" and the monitor will shut off.

Tie the Cuff
1. Plug in the cuff connecting tube into the unit.
2. Wrap the arm cuff, leaving the end of the cuff around the cuff.

3. Insert arm in the cuff. The tube must be aimed in the direction of the hand. Turn the palm of the hand upward and rest the edge of the cuff at about 2-3 cm (0.8\"/>

Helpful Tips
Here are a few helpful tips to help you obtain more accurate readings:
• Do not measure your blood pressure immediately after consuming a large meal.
• To obtain more accurate readings, please wait one hour after meal before measuring.
• Do not smoke or drink tea, coffee or alcohol before measuring your blood pressure.
• Do not measure your blood pressure immediately after taking a bath; wait at least 20 minutes.
• You should not be physically tired or exhausted while taking a measurement.
• Do not measure your blood pressure when you want to discharge urine.
• It is important that you relax during measurement.
• Do not take measurements if you are under stress or tense.
• Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.
• Rest for 5 minutes before measuring.
• Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
• Take the measurement in a silent room.
• Do not move and do not talk when taking a measurement
• Do not cross your legs and keep your feet on the ground.
• Keep your back against the backrest of the chair.
• For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

1. When the monitor is off, press the **START/STOP** to turn on the monitor and start the measurement.
2. Press the **START/STOP** to power off, otherwise it will turn off within 1 minute.
If the **♥** icon appears, please, read carefully the Helpful Tips section and repeat the measurement.
Notes: To stop the measurement at any time, press the START/STOP button under the display. The cuff will quickly deflate. Store the cuff carefully.

Recall the measurements
1. When the Blood pressure monitor is off, select area A or area B and press the "MEM" button to see the average of the last 3 recordings taken within 30 minutes.
If the recorded measurements are less than 3 as well as the last 3 recordings are taken in more than 30 minutes, you will see the latest reading.
2. Press "MEM" or "SET" to see the recording wanted.

WARNING
The most recent recording (1) is shown first. Every time a new recording is made, the first (1) record is assigned to it. All the others are moved by one figure (e.g. 2 becomes 3, and so on), whereas the last record (100) is deleted from the list.

Delete the measurements
If you did not get the correct measurement, you can delete all results by following below steps.
1. When the blood pressure monitor is in memory mode, select area A or area B and keep the "MEM" button pressed for 3 seconds; the display "dE L" blinks. Press "SET" to confirm deleting and the monitor will shut off.
2. If you don't want to delete the records, press "START/STOP" to escape.
3. If there is no record the display will show ":-".

Maintenance
In order to get the best performance from your monitor, please follow the below instructions.
• Put in a dry place and avoid the sunshine.
• Avoid touching water, clean it with a dry cloth in case.
• Avoid the intense shaking and collision.
• Avoid the dusty and unstable temperature environment.
• Using the wet clothing to remove the dirt.
• Avoid washing the cuff.

Calibration and service
• The accuracy of this monitor has been carefully tested and is has been designed for a long life. Indeed the monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or five years of normal use.
• It is generally recommended to check the unit every two years to ensure correct functioning and accuracy. Please consult the Contact Information in this user manual to get in touch with your local authorized Service Center.

WARNINGS
• This device is intended for non-invasive measuring and monitoring of arterial blood pressure on adults, in home and professional environment.
• This device is not suitable for use on newborns, children or pregnant women. The device is not validated to be used on the above listed patient and it could be dangerous. The use on teenagers shall be evaluated by physician.
• The device cannot be used on pre-eclamptic female.
• The device is not suitable for use if your arm circumference is not included in the range 22cm - 42cm (8 3/4\"/>

• The user must check that the equipment functions safely and see that it is in proper working condition before being used.
• Manufacturer will make available on request circuit diagrams, component parts list etc. only to qualified person.
• Too frequent measurements can cause injury to the PATIENT due to blood flow interference. Please check (by observation of the limb concerned) that operation of maxiRAPID does not result in prolonged impairment of the circulation of the blood of the PATIENT.
• This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's arm and fingers will become anaesthetic, swollen and even purple due to a lack of blood.
• Check that the device tubes are not kinked while using the device: kinked tubes can generate lack of blood in the extremities.
• Please use the device under the environment which was provided in the instruction manual. If this is not complied with, correct operation and life of the device may be compromised.
• During using, the patient will contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.
• The application of the CUFF over a wound, as this can cause further injury
• The application of the CUFF and its pressurization on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present because of temporary interference to blood flow and could result in injury to the PATIENT
• The application of the CUFF and its pressurization on the arm on the side of a mastectomy can cause injury or give not accurate measurement
• Do not inflate the cuff on the same limb which other monitoring ME equipment is applied around simultaneously, because this could cause temporary loss of function of the simultaneously-used monitoring ME equipment.
• The device is a non-AP/AVG device and not suitable for use in the presence of anaesthetic mixtures that are inflammable with air, oxygen or nitrous oxide.
• The operator shall not touch DC output jack of AC/DC adapter and the patient simultaneously.
• Please use accessories and detachable parts specified/authorized by the manufacturer. Otherwise, it may cause damage to the unit or danger to the user/patients.
• We recommend checking its performance every two years or after it has been repaired. The device must be checked again in case of damage due to impact (such as falls) or exposure to liquids and/or extreme temperatures (hot/cold) or extreme humidity variations.
• For public place use, please contact your local distributor or the Italian toll-free number 800 900 080 to get support how to proceed.
• Please dispose of accessories, detachable parts, and the ME EQUIPMENT according to the local guidelines.
• In case of problems concerning performances, settings, maintenance or use of the device, contact your local distributor or the Italian toll-free number 800 900 080. Warning! Do not service and/or repair the appliance while it is being used! Any maintenance/service activity/operation must be done solely by a PIC Service Centre. No modification of the equipment are allowed
• When the device is stored at the minimum/maximum storage temperature it takes 4h to warm/cool until the device is ready for its intended use.

This section includes a list of error messages and frequently asked questions for problems you may encounter with your sphygmomanometer. If the products not operating as you think it should, check here before arranging for servicing.

PROBLEM	SYMPTOM	CHECK THIS	REMEDY
No power	Display will not light up.	Batteries are exhausted.	Replace with new batteries
		Batteries are inserted incorrectly. AC adaptor is inserted incorrectly.	Insert the batteries correctly Insert the AC adaptor correctly
Low batteries	The display is dim or shows Lo +	Batteries are exhausted.	Replace with new batteries

Error message	Message	Action
E 1 shows	The cuff is not secure	Refasten the cuff and then measure again.
E 2 shows	The cuff is very tight	Refasten the cuff and then measure again.
E 3 shows	The pressure of the cuff is excess.	Relax for a moment and then measure again.
E10 or E11 shows	The monitor detected motion while measuring.	Movement can affect the measurement. Relax for a moment and then measure again.
E20 shows	The measurement process does not detect the pulse signal.	Loosen the clothing on the arm and then measure again
E21 shows	The treatment of the measurement failed.	Relax for a moment and then measure again.
EExx shows on the display.	A calibration error occurred.	Retake the measurement. If the problem persists, contact the dealer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions.
OUT	Out of range measurement.	Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your Physician.

SPECIFICATIONS

- Power supply:** Battery powered mode: 4 batteries 1.5V AAA; 6V ≈ 315 measurement) Power supply via mains adapter: Input: 100 - 240V~, 50 / 60Hz, 0.2A Max, Output: 6V ≈ 1A Can be powered only with the "Adapter for blood pressure monitor" (REF 02000520000000 (BLJ06L060100P-V)) with European plug and the "Adapter for blood pressure monitor UK" (REF 02010310000000 (BLJ06L060100P-B)) with UK plug (Not included).
- Display:** LCD, 61 x 93mm
- Measurement mode:** Oscillometric inflation -type measurement mode
- Measurement range:** Rated cuff pressure: 0 - 299mmHg
Measurement pressure: SYS 60 - 230mmHg, DIA 40 - 130mmHg
Pulse value: 40 - 199 beat/minute
PRESSURE: ± 3mmHg
Pulse value: ± 5%
- Accuracy:** Temperatures: from 5°C to 40°C. Relative humidity ≤85%
Atmospheric pressure: 86kPa to 106kPa
- Normal working condition:** Temperatures: from -5°C to 50°C. Relative humidity ≤85%
Atmospheric Pressure: 50 - 106 kPa
About 22cm - 42cm
- Storage & transportation condition:** Approx. 250g (Excluding the dry cells)
- Arm circumference:** Approx. 123.5x140x58.5mm
- Net Weight:** Adaptor (REF 02000520000000 with European plug - REF 02010310000000 with UK plug)
- External dimensions:** Continuous operation
- Accessories:** 10,000 Measurements
- Mode of operation:** Type BF applied part (cuff)
- Useful life:**
- Degree of protection:** IP21 degree of protected against solid foreign objects of 12.5mm Ø and greater and protected against vertically falling water drops.
- Protection against ingress of water:** A01
- Software version:**

Contact Information
For further information on our products, please visit the website www.picsolution.com. Here you can find customer service, solve ordinary problems.

EMC Guidance

1) This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
2) "Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
3) Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
4) * Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacturer's declaration – electromagnetic emissions		
The maxiRAPID (REF 02022535000000) is intended for use in the electromagnetic environment specified below. The customer or the user of the maxiRAPID (REF 02022535000000) should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The maxiRAPID (REF 02022535000000) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class [B]	The maxiRAPID (REF 02022535000000) is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The maxiRAPID (REF 02022535000000) is intended for use in the electromagnetic environment specified below. The customer or the user of the maxiRAPID (REF 02022535000000) should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV Input/output lines: ±1 kV	Power supply lines: ±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV line(s) to earth: ±2 kV 100 kHz repetition frequency	line(s) to line(s): ±1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0%U _n ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0%U _n ; 1 cycle and 70%U _n ; 25/30 cycles Single phase: at 0° 0%U _n ; 300 cycle	0%U _n ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0%U _n ; 1 cycle and 70%U _n ; 25/30 cycles Single phase: at 0° 0%U _n ; 300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_n is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity			
The maxiRAPID (REF) 02022535000000 is intended for use in the electromagnetic environment specified below. The customer or the user of the maxiRAPID (REF) 02022535000000 should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	150 kHz to 80 MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150 kHz to 80 MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances: $d = 0.35 \sqrt{P}$ $d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	10V/m, 80% Am at 1kHz	10V/m, 80% Am at 1kHz	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$ 800 MHz to 2,7 GHz $d = 2.3 \sqrt{P}$ Where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
* Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the maxiRAPID (REF) 02022535000000 is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. * Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the maxiRAPID (REF) 02022535000000.			
The maxiRAPID (REF) 02022535000000 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the maxiRAPID (REF) 02022535000000 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=3.5 \sqrt{P}$	80MHz to 800MHz $d=1.2 \sqrt{P}$	800MHz to 2,7GHz $d=2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			


Guidance and manufacturer's declaration - electromagnetic immunity							
The maxiRAPID (REF) 02022535000000 is intended for use in the electromagnetic environment specified below. The customer or the user of the maxiRAPID (REF) 02022535000000, should assure that it is used in such an environment.							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test frequency (MHz)	Band *) (MHz)	Service *)	Modulation *)	Modulation *) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380 – 390	TETRA 400	Pulse modulation *) 18 Hz	1,8	0,3	27
	450	430 – 470	GMRS 460, FRS 460	FM c) ±5 kHz deviation 1 kHz sine	2	0,3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation *) 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation *) 18 Hz	2	0,3	28
	870						
	930						
	1 720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation *) 217 Hz	2	0,3	28
	1 845						
	1 970						
	2 450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation *) 217 Hz	2	0,3	28
	5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation *) 217 Hz	0,2	0,3	9
	5 500						
	5 785						
NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.							

a) For some services, only the uplink frequencies are included.
b) The carrier shall be modulated using a 50% duty cycle square wave signal.
c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = \frac{6}{d} \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

Safety Information			
The symbols below can be found in the Instruction Manual as well as on labels and other components.			
	*PLEASE REFER TO THE INSTRUCTIONS FOR USE* symbol		*TYPE BF APPLIED PART* symbol
	Symbol for "COMPLIES WITH MDD 93/42/EEC as amended by directive 2007/47/EC"		Storage and transport humidity symbol
	*"MANUFACTURER" symbol		Storage and transport temperature symbol
	Symbol for "SERIAL NUMBER"		Caution
	Symbol for "DIRECT CURRENT"		Lot number
	Latex free		Catalogue Number
	Protected against solid foreign objects of 12.5mm Ø and greater and protected against vertically falling water drops.		Atmospheric pressure



Appliance date of production: the second pair of digits in the lot number identifies the year, the third pair identifies the month. (e.g. LOT 1016100001: 16=2016; 10=October).

THE PRODUCT COMPLIES WITH:
EN 1060-1:1995/A2:2009 Non-invasive sphygmomanometers – Part 1: General requirements
EN 1060-3:1997/A2:2009 Non-invasive sphygmomanometers – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems
EN 1060-4:2004 Non-invasive sphygmomanometers – Part 4: Procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers

 **THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.**
The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

 **THIS PRODUCT COMPLIES WITH EU DIRECTIVE 2006/66/EC**
The crossed bin symbol on the batteries indicates that, at the end of their life, they must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for bringing batteries, at the end of their life, to the appropriate collection facilities in order to facilitate treatment and recycling. The correct disposal of your old toy will help prevent potential negative consequences for the environment and human health and favour the recycling of the materials of which the product is made. Illegal disposal of the product by the user will damage the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

Pouch cleaning symbols

-  Wash by hand in cold water,
-  Do not bleach,
-  Do not tumble dry,
-  line drying,
-  Do not iron,
-  Do not dry clean.

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