

OPEN
HERE

TIPS:

This product is disposable.
It is recommended to wear it continuously for 4 hours and then replace it.
Please wear it with the white side facing your face

OPEN
HERE



**NON STERILE
DISPOSABLE MEDICAL FACE MASK**



Warning:

3PLY WITH ELASTIC EAR LOOP


50PCS

1.Failure to follow all instructions and limitations on the use of this product could reduce the efficacy of the respirator and result in sickness or death.
2.A properly selected respirator is essential to protect your health.
Before using this respirator,consult an Industrial Hygienist or Occupational Safety Professional to determine the suitability for your intended use.
3.This product does not supply oxygen.Use only in adequately ventilated areas containing sufficient oxygen to support life.Do not use this respirator when the oxygen concentration is less than 19.5%.

4.Do not use when concentrations of contaminants are immediately dangerous to health or life.
5.Leave work area immediately and return to fresh air if:(a) breathing difficult (b) dizziness or other distress occurs.
6.Facial hair,beards,and certain facial Characteristics may reduce the effectiveness of this respirator.
7.Never alter or modify this respirator in any way(except as indicated in the instructions).
8.This product is intended for single use only.No maintenance necessary.Discard used respirator after single use.
9.Keep respirators in the display box away from direct sunlight until use.


Color: white blue green pink yellow

How to wear:
white side of the mask
facing inside,color side
of the mask facing outside



Executive Standard:EN14683:2019
BFE≥95 TYPE I

50PCS



3PLY WITH ELASTIC EAR LOOP



**NON STERILE
DISPOSABLE MEDICAL FACE MASK**



**NON STERILE
DISPOSABLE MEDICAL FACE MASK**



3PLY WITH ELASTIC EAR LOOP

Directions for use:

Pull the elastic ear loop around the ear and position for best fit.

Special advantages:

- Extremely lightweight.
- Does not impair breathing,voice and vision.
- No interference with glasses.



EC REP

Hospital & Homecare Medical Device Co., Ltd. China.
No. 12 Plant, Phase II, Standard Plant,
Gaochun Economic Development Zone,Nanjing

CMC Medical Devices & Drugs S.L
C/Horacio Lengo Nº 18
CP 29006, Málaga-Spain
Tel: +34951214054
Fax: +34952330100
email - info@cmcmedicaldevices.com



**NON STERILE
DISPOSABLE MEDICAL FACE MASK
3PLY WITH ELASTIC EAR LOOP**

- HYPOALLERGENIC
- VERY LOW RESISTANCE TO BREATHING
- NOSE BAR ADAPTABLE
- HIGH FILTRATION CAPACITY
- PERFECT FITTING

UDI-DI:

